

## Poly-HRP IHC Amplification Kit

### Introduction

CytoMol's Poly-HRP IHC Amplification Kit utilizes a novel controlled polymerization technology to prepare polymeric HRP-linker antibody conjugates. Comparing to conventional biotin-streptavidin based detection kits; the Poly-HRP IHC Amplification Kit has the advantages of simpler protocols, high amplification power, biotin free and more consistent immunostaining outcomes on archival tissues and on difficult-to-work antibodies. These advantages would bring to laboratories the benefit of more accurate result, faster turn-around, less trouble-shooting and better cost saving.

### Features

- High sensitivity
- Simple protocol

### Application

- Immunohistochemistry

### Description

Affinity purified anti mouse/rabbit IgG, polymeric HRP linked antibody is enclosed in Poly-HRP IHC Amplification Kit. With the special ready-to-use Blocking Solution and 30x DAB Solution A and B, this kit generates high quality immunohistochemistry results.

### Quality Control

In order to guarantee the kits' quality, each lot of the kit is tested with primary antibody to show the positive staining on control slides.

### Contents

Item	Size	Storage
Poly-HRP Anti Mouse/rabbit IgG	15 ml	2-8°C
Blocking Solution	15 ml	2-8°C
DAB Solution A 30x	1 ml	2-8°C
DAB Solution B 30x	1 ml	2-8°C

### IHC Protocol: (at room temperature)

Deparaffinization, rehydration

Blocking endogenous peroxidase with 3% H<sub>2</sub>O<sub>2</sub> in dd water for 10 minutes

Performing antigen retrieval or enzyme digestion if necessary

1. Block with Blocking Solution\*, 10 min; Blot off blocking solution but no need to wash
2. Incubate with mouse or rabbit primary antibody, 30-60 min, wash with TBS, 5'x2
3. Incubate with Poly-HRP anti-Mouse/Rabbit IgG, 30 min, wash with TBS, 5'x2
4. Incubate with DAB solution\*\*, 5 min, wash with dd water
5. Counter stain with hematoxylin and mount

#### Note:

\* The blocking solution can also be used as diluent for primary antibodies. When blocking solution is used as the diluent for primary antibody, blocking step can be omitted.

\*\* Instruction: Mix one drop DAB Solution A (33 µl) and one drop DAB Solution B (33 µl) with 0.93 ml dd water, then apply to tissue sections